

# JAAS

Accepted Manuscript



This is an *Accepted Manuscript*, which has been through the Royal Society of Chemistry peer review process and has been accepted for publication.

*Accepted Manuscripts* are published online shortly after acceptance, before technical editing, formatting and proof reading. Using this free service, authors can make their results available to the community, in citable form, before we publish the edited article. We will replace this *Accepted Manuscript* with the edited and formatted *Advance Article* as soon as it is available.

You can find more information about *Accepted Manuscripts* in the [Information for Authors](#).

Please note that technical editing may introduce minor changes to the text and/or graphics, which may alter content. The journal's standard [Terms & Conditions](#) and the [Ethical guidelines](#) still apply. In no event shall the Royal Society of Chemistry be held responsible for any errors or omissions in this *Accepted Manuscript* or any consequences arising from the use of any information it contains.

1  
2  
3 **Direct mercury analysis in environmental solids by ICP-MS with**  
4 **on-line sample ashing and mercury pre-concentration using a direct**  
5  
6  
7  
8  
9  
10 **mercury analyzer**  
11  
12  
13  
14  
15  
16

17 Derek D. Bussan, Ryan F. Sessums, James V. Cizdziel\*  
18  
19  
20  
21  
22  
23

24 Department of Chemistry and Biochemistry, University of Mississippi, University, Mississippi  
25  
26 38677, USA  
27  
28

29 \*Corresponding author; email: [Cizdziel@olemiss.edu](mailto:Cizdziel@olemiss.edu); Tel: 662-915-1814  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

**Abstract**

A commercially available direct mercury analyzer (DMA) based on sample combustion, pre-concentration by amalgamation with gold, and atomic absorption spectrophotometry (AAS) was coupled to a sector field inductively coupled plasma mass spectrometer (ICP-MS). The combined system (DMA-ICP-MS) was optimized and evaluated for the determination of mercury in environmental samples using both external calibration and isotope dilution for quantitation. The method was validated using certified reference materials, including sediment, leaves, and fish-muscle tissue. The limit of detection was ~0.37 pg, about two orders of magnitude lower than the DMA system alone. Precision was generally <7% relative standard deviation, and total analyses time per sample was <8 min. The DMA-ICP-MS system has several advantages over both the DMA alone and conventional cold-vapor AAS for the determination of mercury, including increased sensitivity, lower detection limits, decreased potential for sample contamination, and applicability to Hg stable isotope tracer studies.

**Keywords:** mercury, direct mercury analysis, ICP-MS, Hg isotopes, thermal decomposition, isotope dilution

## Introduction

Mercury (Hg) is a persistent, mobile, and highly toxic environmental pollutant that is routinely measured in environmental and toxicological studies and biomonitoring programs.<sup>1-3</sup> Direct mercury analyzers (DMAs) are commonly used for determining total-Hg in a variety of sample matrices. These instruments integrate sample combustion, matrix removal and Hg pre-concentration via amalgamation with gold, and atomic absorption spectrophotometry. They are popular because they increase analytical throughput, reduce waste generation and the potential for sample contamination, and have relatively low operating and capital costs. However, atomic absorption spectrometry (AAS) and atomic fluorescence spectrometry (AFS) based methods are unable to distinguish individual isotopes, and thus are not suitable for quantitation using isotope dilution or for isotope tracing studies.

Isotope dilution mass spectrometry (IDMS) is an accurate method to quantify elements that have more than one isotope, and it has long been used for the determination of Hg in environmental samples.<sup>4</sup> In IDMS, a sample is mixed with an isotopic standard enriched in one of the minor isotopes of the analyte element, and, after equilibration, mass spectrometry is used to measure the altered isotope ratios. The resulting isotopic compositions, along with the mass of the sample and spike, are used to calculate the concentration of the element in the sample. Mercury stable isotopes have also been used as tracers to investigate Hg transformations in the environment, and, for more than a decade, these studies have yielded significant new insight into the processes controlling Hg speciation and bioavailability.<sup>5</sup>

Coupling a DMA to an inductively coupled plasma mass spectrometer (ICP-MS) offers several potential advantages over both the DMA alone and conventional cold vapor (CV)-AAS or CV-AFS for the determination of Hg, including: (1) increased sensitivity and lower detection

1  
2  
3 limits, (2) efficient (ca. 100%) Hg introduction compared to conventional liquid nebulization, (3)  
4 simple interface design and operation, (4) improved accuracy through use of isotope dilution, (5)  
5  
6 increased speed and throughput of analyses, and (6) the possibility of direct analysis of solid  
7  
8 samples as part of Hg stable isotope tracer studies of environmental processes.  
9  
10

11  
12 Because there are no reports in the literature specifically evaluating a DMA-ICP-MS  
13 coupled system, the aim of the current study was to: (1) develop a simple procedure to couple the  
14 DMA with an ICP-MS, (2) optimize the combined system, (3) evaluate common analytical  
15 figures-of-merit using both external calibration and IDMS, and (4) apply the method to  
16 environmental and biological sample matrices.  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26

## 27 **Materials and methods**

### 28 *Standards and reference materials*

29  
30  
31  $^{201}\text{Hg}$  (93.08%) enriched stable isotope was obtained from Oak Ridge National  
32 Laboratory as HgO. The oxide was dissolved in 10% high-purity (Optima grade)  $\text{HNO}_3$  and  
33 diluted. A single element natural isotopic abundance Hg standard obtained from Spex Certiprep  
34 (Metuchen, NJ, USA) was used to optimize the ICP-MS prior to analysis, and to calculate the  
35 concentration of the enriched stable isotope spike using reverse-IDMS. For the latter, approx.  
36 1.05 g of the natural isotopic abundant Hg (concentration  $493.9 \text{ ng g}^{-1}$ ) was accurately weighed  
37 and added to 0.25 g of the  $415 \text{ ng g}^{-1}$  spike, and the mixture diluted to  $\sim 50.00 \text{ g}$  in 10%  $\text{HNO}_3$ .  
38 The concentration of the  $^{201}\text{Hg}$  spike was determined by reverse IDMS using the same equations  
39 for IDMS.<sup>6</sup>  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51

52  
53 DORM-3 (dogfish muscle), NIST 2709 (sediment), and NIST 1547 (peach leaves) was  
54 used as reference materials. DORM-3 was obtained from the National Research Council of  
55  
56  
57  
58  
59  
60

1  
2  
3 Canada (Ottawa), and the other reference materials were from the National Institute of Standards  
4 and Technology (Gaithersburg, MD, USA). Ultrapure water ( $>18 \text{ M}\Omega \text{ cm}^{-1}$ ) was obtained from  
5  
6 a Barnstead Nanopure Diamond system.  
7  
8

### 9 10 *Direct Mercury Analyzer*

11  
12 A DMA-80 (Milestone Inc., Shelton, CT) was used throughout this study. DMAs and  
13 their applications in environmental and biological studies have been described elsewhere.<sup>7-9</sup>  
14  
15 Briefly, samples were weighed in nickel or quartz boats and placed in an autosampler which  
16 sequentially inserts them into the DMA's combustion tube where they were heated to  $\sim 650 \text{ }^\circ\text{C}$   
17 with oxygen as a carrier gas. Gaseous combustion products were carried through a heated  
18 catalyst, where Hg species are converted to elemental Hg vapor ( $\text{Hg}^0$ ), and where halogens and  
19 other species which can interfere with the analysis are trapped. Mercury vapor and other  
20 decomposition products were carried to a glass tube gold-coated sand where  $\text{Hg}^0$  is selectively  
21 trapped. Later, the trap is rapidly heated to release  $\text{Hg}^0$  vapor into a single beam  
22 spectrophotometer. Mercury concentration is calculated based on the absorbance measured at  
23 253.7 nm and the weight of the sample. In the current study, operating times for drying,  
24 combustion, and post-combustion flushing periods were 60, 180, and 45 s, respectively, for a  
25 total analysis time of  $<5$  min per sample. DMA instrumental settings are given in Table 1.  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42

### 43 *ICP-MS*

44  
45 A sector field-ICP-MS (Element-XR, Thermo Scientific) was used to measure Hg  
46 isotopes. The ICP-MS was optimized prior to analysis using a Hg tuning solution (Spex  
47 Certiprep; Metuchen, NJ, USA) introduced using a microflow concentric nebulizer and a peltier-  
48 cooled cyclonic spray chamber ( $\text{PC}^3$ ; Elemental Scientific, Omaha, NE, USA). Argon from the  
49 DMA was introduced downstream of the  $\text{PC}^3$  just prior to the torch using a T-junction (Figure 1).  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 This setup allows for easy instrument tuning and, if desired, introduction of solutions for  
4 measuring mass discrimination. When tuning the ICP-MS, argon from the DMA was kept on.  
5  
6 When running samples on the DMA-ICP-MS, the optimized nebulizer gas flow was kept on but  
7  
8 the uptake line from the solution was removed.  
9  
10

11  
12 ICP-MS instrumental settings are given in Table 1. Mercury isotopes ( $^{201}\text{Hg}$  and  $^{202}\text{Hg}$ )  
13 were monitored in low resolution mode. A higher resolution setting was not necessary because  
14  
15 interferences are not present as the Hg is isolated from the matrix via the amalgamation step.  
16  
17 The mass window was set to 5% with 100 points per peak, yielding 5 data points (averaged per  
18  
19 pass) across the middle portion of the flat-topped peak. Fast electric scanning (123 runs/pass)  
20  
21 was used with an integration time of 10 ms, resulting in a data point for each isotope every 0.86 s  
22  
23 over a total acquisition time of 1.72 min.  
24  
25  
26  
27  
28

### 29 *DMA combination with ICP-MS: interface and operation*

30  
31 A schematic of the DMA-ICP-MS interface showing the operation of the system during a  
32  
33 typical analytical run is shown in Figure 1. Instrumental settings for both the DMA and the ICP-  
34  
35 MS are given in Table 1. Valves “upstream” and “downstream” of the DMA were used to direct  
36  
37 gases compatible with DMA and ICP-MS. The operation of the valves and triggering of the  
38  
39 ICP-MS was done manually; however, this could certainly be automated by researchers or  
40  
41 instrument manufacturers. To interface the DMA with the ICP-MS, a 3.175 mm (OD) teflon  
42  
43 tube was inserted directly into the outlet of the DMA analysis cell making a tight seal without  
44  
45 any further manipulation. The other end of the tube was connected to a 4-port stainless steel  
46  
47 valve (Valco Inc., Houston, TX, USA). The entire length of the tube was wrapped in silicone  
48  
49 rubber-encapsulated heat tape (Thermolyne-Barnstead) to prevent condensation in the line. The  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 “downstream” valve was used to direct gases leaving the DMA either to vent or to the ICP-MS  
4  
5 via another 3.175 mm (OD) Teflon tube.  
6  
7

8         At the start of a run and during the sample drying and decomposition (sample ashing and  
9 Hg pre-concentration) phase of the analysis, the upstream valve was positioned to deliver O<sub>2</sub> to  
10 the DMA, while the downstream valve was positioned to direct combustion products to vent.  
11  
12 The valves were switched at the conclusion of the decomposition phase (during the start of purge  
13  
14 1 phase, before the amalgamator is heated) to change the carrier gas from O<sub>2</sub> to argon and to  
15  
16 change from vent to the ICP-MS. Doing this early in this stage allows the system time to  
17  
18 stabilize. Next, data acquisition was initiated for the ICP-MS about 10 sec before the  
19  
20 amalgamator was heated to drive off Hg. This analytical scheme allows Hg (free of combustion  
21  
22 products) to be carried to the plasma in argon.  
23  
24  
25  
26  
27  
28

### 29 *DMA-ICP-MS (external calibration)*

30  
31         For external calibration, a single isotope (<sup>202</sup>Hg) was monitored in low resolution mode  
32 and peak areas were used for quantitation. The system was calibrated using a temperature-  
33 controlled Hg-vapor calibration unit (Tekran 2505) and standard solutions. For calibration in the  
34 low range, a gas-tight syringe (Hamilton Company) was used to deliver approximately 50, 100,  
35  
36 150, 200, 250, 300 and 350 pg of Hg to the DMA via an injection port that was outfitted to the  
37  
38 carrier gas stream just prior to entering the combustion chamber (Figure 1). For calibration at  
39  
40 higher levels, a standard solution was weighed into quartz boats using an analytical balance with  
41  
42 0.1 mg readability, and the boats delivered to the DMA. The standard solutions were prepared  
43  
44 from a 10,000 mg L<sup>-1</sup> stock solution (SPEX CertiPrep) and delivered approximately 1.5, 4.0 and  
45  
46 10 ng of Hg to the DMA. The weights of the solutions introduced to the DMA ranged from  
47  
48 about 20 to 50 mg.  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

### *DMA-ICP-MS (isotope dilution quantitation)*

For isotope dilution measurements (single spike method), a  $^{201}\text{Hg}$  enriched isotope standard solution (described above) was used. The samples were accurately weighed in a combustion boat using the same balance described above. The spike solution was then added directly on the sample and the mixture was re-weighed. Spikes were added to deliver  $\sim 4.65$  ng  $^{201}\text{Hg}$  for NIST 2709 and DORM-3, and 0.11 ng for NIST 1547. The boats were delivered to the DMA and the DMA run as described before.  $^{201}\text{Hg}$  and  $^{202}\text{Hg}$  were monitored by the ICP-MS in low resolution mode and peak areas were used to determine the Hg isotope ratio in the sample. The ratio, together with the mass of the sample and spike, was used to quantify total-Hg in certified reference materials and lake sediment using standard isotope dilution calculations.<sup>4,6</sup>

### **Results and discussion**

For conventional IDMS of solution digests, an intimate mix of the isotope spike with the sample is critical to achieve accurate data. This is not as important for samples that are subject to thermal decomposition as long as all the Hg is quantitatively transferred and measured. Moreover, because Hg isotopes will fractionate as they desorb from a heated gold trap (lighter isotopes preferentially desorbing first), it is important that the entire area of the peak be used in the calculations so that the data captures accurate isotope ratios.<sup>10</sup> A typical intensity versus time chromatogram generated by the DMA-ICP-MS is shown in Figure 2. There is very little time delay ( $< 2$  s) from when the DMA amalgamator is heated to when the ICP-MS Hg signal rises.

### *Limit of detection, accuracy, and precision*

For the DMA-alone, the limit of detection (LOD,  $3\sigma$  criteria) was 13 pg based on eight blank measurements (empty sample boats) and the slope of the calibration curve in the low

1  
2  
3 concentration range (gas phase standard injections). Using the same approach the LOD for the  
4  
5 DMA-ICP-MS was 0.37 pg, nearly two orders of magnitude lower than the DMA alone. The  
6  
7 DMA-ICP-MS LOD corresponds to 1.2 pg g<sup>-1</sup> for 250 mg of sediment (a typical weight for  
8  
9 sediment in a DMA boat) or 5.0 pg g<sup>-1</sup> for 60 mg of biological tissue (typical weight for fish  
10  
11 tissue in a DMA boat). Thus, the DMA-ICP-MS method is capable of ultra-trace level analyses.  
12  
13 Moreover, there was virtually no carryover between samples for sample weights that deliver less  
14  
15 about 40 ng of Hg. Nevertheless, because Hg can adsorb on Teflon surfaces, carryover should  
16  
17 be monitored, particularly following samples containing relatively high levels of Hg.  
18  
19  
20  
21

22 The method was validated using a variety of certified reference materials, including  
23  
24 sediment (NIST 2709), leaves (NIST 1547), and fish tissue (DORM 3). Recoveries for these  
25  
26 certified reference materials were within the certified range (Table 2). Following U.S.  
27  
28 Environmental Protection Agency Method 7473, the DMA alone yielded recoveries ranging  
29  
30 from 100-103% of certified values, indicating that the DMA was operating properly. The  
31  
32 corresponding DMA-ICP-MS data was also in good agreement with certified values: ±8% for  
33  
34 external calibration and slightly better accuracy (±2%) for isotope dilution measurements.  
35  
36  
37  
38

39 For external calibration of the DMA-ICP-MS, precision was <7% relative standard  
40  
41 deviation (RSD) (Table 2). For IDMS quantitation, the precision was <4% RSD for both  
42  
43 DORM-3 (n=9) and NIST 2709 (n=5), but higher for NIST 1547 leaves (11% RSD), which had  
44  
45 the lowest Hg levels of the three (31 ± 7 ng g<sup>-1</sup>). Uncertainty in the results can stem from the  
46  
47 relatively small weights used for either the spike solutions or the sample, where a small  
48  
49 difference in weight can have a big impact on results. We kept the absolute amounts of Hg to  
50  
51 <40 ng (~0.5 – 40 ng) to minimize carryover and prevent overloading the DMA with Hg. Thus,  
52  
53 sample weights for NIST 2709 soil (1400 ± 80 ng g<sup>-1</sup>), DORM-3 (382 ± 60 ng g<sup>-1</sup>) and NIST  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 1547 ( $31 \pm 7 \text{ ng g}^{-1}$ ) were ~25 mg, ~125 mg, and ~200 mg, respectively. Also, whereas the  
4  
5 DMA itself and the transfer line from the DMA to the downstream valve was heated to prevent  
6  
7 condensation, we still wanted to keep the volume of liquid (mass of spike) to a minimum (spike  
8  
9 levels were 20-50 mg). Using weights lower than these for the sample and spike are likely to  
10  
11 increase uncertainty in the analytical results.  
12  
13

14  
15 We are currently using the DMA-ICP-MS method to determine total-Hg in sediments  
16  
17 from an old-growth cypress wetland at Sky Lake, an oxbow lake in the Mississippi Delta alluvial  
18  
19 flood plain, as part of a stable isotope tracer study examining Hg methylation/demethylation  
20  
21 rates. So far, the RSD has been <10% (2.8 – 12.6%), based on seven different sediment samples  
22  
23 run in triplicate and ranging in concentration from 35 – 147  $\text{ng g}^{-1}$  (Table 3).  
24  
25  
26  
27  
28

## 29 **Conclusions**

30  
31 A direct mercury analyzer was coupled with an ICP-MS for the direct determination of  
32  
33 Hg in environmental and biological samples. The approach eliminates the need for acid-  
34  
35 digestion of samples and allows for isotope ratio measurements. Compared to both the DMA  
36  
37 alone and conventional CV-AAS, the combined system was shown to lower detection limits,  
38  
39 increase sensitivity, and provide accurate quantification by the isotope dilution method.  
40  
41 Anticipated environmental applications of the method include evaluating Hg burdens in  
42  
43 suspended sediments, non-lethal fish-Hg determinations using milligram quantities of tissue,  
44  
45 exposure studies using temporally-resolved hair analysis, and environmental process studies  
46  
47 using Hg stable isotope tracers.  
48  
49  
50  
51  
52  
53  
54

## 55 **Acknowledgements**

1  
2  
3 We thank Dr. Edward Heithmar and Ms. Charlita Rosal (U.S. Environmental Protection  
4 Agency) who, along with the corresponding author, presented the analytical approach behind the  
5 present work at Winter Conference on Plasma Spectrochemistry in 2002. We also thank Oscar  
6 “Beau” Black, Kristen Dickerson, Ruiqi Feng, and Jingjing Chen for help with various aspects of  
7 the project. The ICP-MS used in this study was obtained through a U.S. National Science  
8 Foundation grant (Award #0923080).  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

**Table 1. DMA-ICP-MS Instrumental Settings**

DMA Parameters	
Gas flow	200 mL min <sup>-1</sup>
Drying	200 °C for 60 s
Decomposition	650 °C for 180 s
Purge 1	60 s
Amalgamator heat	~850 °C for 12 s
Purge 2	45 s

Plasma Parameters	
Cool gas flow	16 L min <sup>-1</sup>
Aux. gas flow	1.03 L min <sup>-1</sup>
Sample gas flow	1.05 L min <sup>-1</sup>
RF power	1280 W

ICP-MS Data Acquisition Parameters	
Isotopes	<sup>201</sup> Hg, <sup>202</sup> Hg
Resolution	Low
Mass Window	5%
Points per peak	100

**Table 2. Total-Hg (ng g<sup>-1</sup>) in reference materials determined by thermal decomposition, amalgamation, atomic absorption spectrometry (Method 7473) and by the DMA-ICP-MS using external calibration or isotope dilution mass spectrometry**

Matrix	Reference Material	Certified	Observed								
			Method 7473			External Calibration			IDMS		
			Conc.	RSD (%)	Recovery (%)	Conc.	RSD (%)	Recovery (%)	Conc.	RSD (%)	Recovery (%)
Sediment	NIST 2709	1400 ± 80	1396 ± 23	1.6	100	1510 ± 80	5.3	108	1416 ± 24	1.7	101
Leaves	NIST 1547	31 ± 7	32 ± 2	5.2	103	29 ± 2	6.9	94	31 ± 4	11	100
Fish-muscle	DORM-3	382 ± 60	387 ± 23	5.9	101	420 ± 15	3.6	107	373 ± 12	3.3	98


n = 5, except for DORM-3 where n = 9

**Table 3. Total-Hg in sediment from Sky Lake determined by DMA-ICP-MS using IDMS**


Sample	Date	T-Hg (ng g <sup>-1</sup> ) (n=3)	RSD (%)
Wetland	11/30/2014	88.7 ± 4.4	5.0
Wetland		91.9 ± 6.4	6.9
Open Water		35.5 ± 3.2	9.0
Wetland	1/25/2015	147 ± 4	2.8
Wetland		119 ± 15	12.6
Open Water		110 ± 5	5.0
Open Water		106 ± 8	7.7

## References


1. O. Malm, *Environmental Research*, 1998, **77**, 73-78.
2. P. B. Tchounwou, W. K. Ayensu, N. Ninashvili and D. Sutton, *Environmental Toxicology*, 2003, **18**, 149-175.
3. R. F. Palmer, S. Blanchard, Z. Stein, D. Mandell and C. Miller, *Health & Place*, 2006, **12**, 203-209.
4. R. G. Smith, *Analytical Chemistry*, 1993, **65**, 2485-2488.
5. H. Hintelmann, K. Keppel-Jones and R. D. Evans, *Environmental toxicology and chemistry*, 2000, **19**, 2204-2211.
6. J. I. G. Alonso and P. Rodriguez-Gonzalez, *Isotope dilution mass spectrometry*, Royal Society of Chemistry, 2012.
7. C. T. Costley, K. F. Mossop, J. R. Dean, L. M. Garden, J. Marshall and J. Carroll, *Analytica Chimica Acta*, 2000, **405**, 179-183.
8. J. V. Cizdziel, T. A. Hinnens and E. M. Heithmar, *Water, Air, and Soil Pollution*, 2002, **135**, 355-370.
9. EPA, Washington, DC, 1998, pp. SW846, Method 7473, <http://www.epa.gov/osw/hazard/testmethods/sw7846/pdfs/7473.pdf>.
10. R.D. Evans, H. Hintelmann and P.J. Dillon, *Journal of Analytical Atomic Spectrometry*, 2001, **16**, 1064-1069.



430x202mm (96 x 96 DPI)



Schematic of the DMA-ICP-MS interface shown in ICP-MS analysis mode (left) and operation of the valves during a typical run (right).  
430x202mm (96 x 96 DPI)



Typical intensity vs. time chromatogram generated by the DMA-ICP-MS. At ~10 sec in the figure the DMA gold trap is heated releasing Hg to the ICP-MS. Numbers in the figure represent four of the stable isotopes of Hg that were monitored.  
197x136mm (96 x 96 DPI)