
Volume 1 Number 4 October 2020 Pages 159-290

RSC Chemical Biology



ISSN 2633-0679



COMMUNICATION

Kousuke Tsuchiya, Yutaka Kodama, Keiji Numata *et al.* Cellular internalization mechanism of novel Raman probes designed for plant cells

RSC Chemical Biology



COMMUNICATION

View Article Online
View Journal | View Issue



Cite this: *RSC Chem. Biol.*, 2020, **1**, 204

Received 17th July 2020, Accepted 29th July 2020

DOI: 10.1039/d0cb00128g

rsc.li/rsc-chembio

Cellular internalization mechanism of novel Raman probes designed for plant cells†

Keiko Midorikawa,‡^a Kousuke Tsuchiya, \$\bar{\mathbb{D}}\disp**\text{ab}\$ Simon Sau Yin Law,\$^a Yu Miyagi,\$^a Takuya Asai,\$^c Takanori lino, \$\bar{\mathbb{D}}^c\$ Yasuyuki Ozeki, \$\bar{\mathbb{D}}^c\$ Yutaka Kodama \$\bar{\mathbb{D}}^**\text{ad}\$ and Keiji Numata \$\bar{\mathbb{D}}^**\text{ab}\$

Diphenylacetylene derivatives containing different polymeric components, poly(L-lysine) (pLys) or tetra(ethylene glycol) (TEG) were designed as novel Raman imaging probes with high Raman sensitivity and low cytotoxicity in living plant cells. The pLys-conjugated probe is internalized *via* an endocytosis-dependent pathway, whereas TEG-conjugated probe most likely induces direct penetration into the plant cells.

Plant cells synthesize many more secondary metabolites than animal cells, many of which are useful to humans. 1,2 Therefore, much effort has been made to stabilize production and increase vields for these metabolites by plant biotechnology. Recently, an intracellular delivery system using a complex with a polymer carrier has emerged as a tool that enables cellular internalization of various bioactive molecules.^{3,4} These molecules contribute to the improvement of production efficiency of metabolite with genetic modification. In order to further improve the delivery of materials it is important to elucidate the internalization mechanisms that remain unknown. Thus, the tools are needed to visualize the process of target molecule internalization in living cells. Site-selective staining using fluorescent probes has been developed to visualize how biomolecules are localized in living cells using fluorescence microscopy. 5 However, following the target molecules themselves is difficult in fluorescence microscopy using conventional probe dyes unless the target molecules are chemically labelled with a fluorescent

Notably, the metabolites of interest must be independently modified with a Raman-detectable moiety for visualization in the previous examples. To date, analysis using Raman microscopy has been applied effectively to both animal and plant cells. 9-12 Therefore, a novel Raman-based monitoring system would be highly beneficial for the visualization of various phenomena associated with molecular translocation in living plant cells. However, in general, Raman spectromicroscopic analysis has lower detection sensitivity than fluorescence probe imaging. To enhance the sensitivity, Raman probes require a chemical structure with a stronger Raman signal intensity, such as diphenylacetylene moieties. The drawbacks of these structures are hydrophobic, water-insoluble, cytotoxic, and poorly permeable into plant cells. 13,14 Hence, we have overcome these problems by adding a hydrophilic and highly biocompatible polymeric carrier to diphenylacetylene, which shows strong Raman intensity. Besides, further functionalization enables such the polymeric component-fused Raman probes to be utilized for delivery of various target molecules into a cell, as we have reported that the polymeric carriers can be easily

probe via in vivo reactions. On the other hand, Raman microscopy can detect a vibration attributed to a specific chemical structure within target molecules and can be applied to observe molecular behaviour in living cells without invasion or staining. The alkyne group is the best candidate as a Raman probe for modifying target molecules because it produces a specific peak in a silent region of the Raman spectrum, 6 which enables detection of the probe in molecular crowding situations found in living cells. Furthermore, since the acetylene moiety is a bioorthogonal functional group with an inert nature to chemical reactions in vivo, it has the advantage that it is less susceptible to non-specific interactions and intracellular reactions. Therefore, the fluorescent probe needs to be designed in consideration of non-specific interactions in vivo and unexpected chemical reactions, whereas the Raman probe can be applied to a wide range of targets. In recent years, the introduction of alkynyl groups into low-molecular-weight compounds, such as drugs and metabolites, has been reported for visualization in vivo.7,8

^a Biomacromolecules Research Team, RIKEN Center for Sustainable Resource Science, 2-1 Hirosawa, Wako, Saitama 351-0198, Japan. E-mail: keiji.numata@riken.jp, kosuke.tsuchiya@riken.jp

b Department of Material Chemistry, Graduate School of Engineering, Kyoto University, Kyoto-Daigaku-Katsura, Nishikyo-ku, Kyoto 615-8510, Japan

^c Department of Electrical Engineering and Information Systems, The University of Tokyo, Tokyo 113-8656, Japan

^d Center for Bioscience Research and Education, Utsunomiya University, Tochigi 321-5805, Japan. E-mail: kodama@cc.utsunomiya-u.ac.jp

[†] Electronic supplementary information (ESI) available. See DOI: 10.1039/d0cb00128g

[‡] These authors are co-first author.

combined with arbitrary functional molecules via the thiol-ene reaction when the polymeric carriers are modified with a maleimide group.3

Herein, we developed novel Raman probes that enable us to monitor the active translocation of target molecules in plant cells with high Raman sensitivity and low cytotoxicity to plant cells. We combined diphenylacetylene with two types of polymeric components, poly(L-lysine) (pLys) and poly(ethylene glycol) (PEG) of different lengths. The pLys and PEG components, which are both hydrophilic but membrane-interactive or physiologically inert, respectively, are often used as part of carriers to deliver exogenous molecules. The two types of conjugated materials were found to function as efficient Raman probes in plant cells with different uptake mechanisms depending on the polymeric components. Furthermore, to improve the sensitivity, we used stimulated Raman scattering (SRS) microscopy. 15,16 The current system combining the novel Raman probes and SRS imaging allows us to analyse molecular behaviours in live cells at a practical level.

Amine-terminated pLys (n = 4-8, the average DP is 6), tetra(ethylene glycol) (TEG), and PEG (n = 24) were conjugated with 4-phenylethynylphthalic anhydride to obtain the Raman probes, PhC≡CPh-pLys (1), PhC≡CPh-TEG (2), and PhC≡CPh-PEG (3), respectively (Fig. 1a). The peptide pLys was selected as membrane-permeable component according to the screening study of cell penetrating peptides (CPPs) for plant cells. ¹⁷ The CPPs designated for animal cells showed different internalization behaviour in plant cells, and Lys-containing peptides showed a relatively high internalization efficiency in plant cells.

The chemical structures of the probes were confirmed by ¹H NMR spectroscopy and matrix-assisted laser time-of-flight (MALDI-TOF) mass spectrometry (Fig. S1-S6, ESI†). The alkynetagged thymidine analog 5-ethynyl-2'-deoxyuridine (EdU) was used as an internal standard for the probe solutions, 18 and the obtained spectra were normalized to compare the Raman intensities. All three probes showed a strong peak in the same Raman shift region (2215 cm⁻¹). Additionally, we confirmed that the intensities of the three probes were similar (Fig. 1b), with higher intensity than the standard ethynyl compound EdU.

Although both EdU and diphenylacetylene have alkynyl groups, diphenylacetylene shows higher Raman intensity because of the aromatic rings at both ends, which leads to the enhancement of Raman scattering. We confirmed that this enhanced Raman effect was maintained when the polymeric component was combined. Then, in order to investigate whether these probes are applicable to plant cells, we examined the internalization of the probes by Raman microscopy using tobacco bright yellow 2 (BY-2) cells and Arabidopsis thaliana. The BY-2 cell line is the most widely used plant cell culture. 19 After the BY-2 cells were incubated in a probe-containing medium, the peak at 2215 cm⁻¹ assignable to the diphenylacetylene moiety was detected within the cytoplasm of cells in the Raman spectra (Fig. 2a and b).

The internalization of 1 and 2 into BY-2 cells was observed, whereas no Raman peak for the diphenylacetylene moiety was detected in the cells treated with 3. We confirmed that the




Fig. 1 Relative Raman intensity of new probes. (a) Chemical structure of Raman probes. PhC≡CPh-pLys (1), PhC≡CPh-TEG (2), PhC≡CPh-PEG (3). (b) Comparison of Raman intensity/shift between three types of probes. The spectrum was normalized to the Raman intensity of EdU.

signal from Raman probe 1 and 2 was detected at multiple points in one BY-2 cell, mainly within the cytoplasm region (Fig. S7 and S8, ESI†). These results indicate that the pLys and TEG components function as effective carriers to internalize the hydrophobic diphenylacetylene moiety, which hardly dissolves in an aqueous buffer. In contrast, the long PEG chain hampered the internalization of the probe, probably because a long PEG chain shows minimal nonspecific interaction with biomolecules, including membrane lipids, due to its large exclusion volume.20 Similarly, we carried out an experiment using Arabidopsis thaliana, which is a model plant. In the experiment with leaf cells, we used var2-1 mutant plant having variegated leaves to suppress autofluorescence derived from chlorophyll in leaves (Fig. S9, ESI†). 21-23 The var2-1 mutant plant has mutations in the subunits of the thylakoid-localized FtsH complex that are likely involved in chloroplast biosynthesis and development. 24,25 The leaves were dipped in the 100 μM probe solution, then washed and observed with a Raman microscope. As a result, internalization into cytoplasm was confirmed for 1





Fig. 2 Internalization of Raman probes into BY-2 cells and Arabidopsis thaliana leaf cells. (a) Squares indicate irradiated points in the cytoplasm in BY-2 cells (1 μ m imes 1 μ m in size). Scale bars are 10 μ m. (b) Averaged Raman spectra in the cytoplasm of BY-2 cells were obtained after 5 h treatment with 100 μ M **1**, **2**, or **3**. Diphenylacetylene shows a Raman peak at 2215 cm⁻¹. (c) Squares indicate irradiated points in the cytoplasm (cp) and vacuole (vc) in a var2-1 leaf cells (1 μ m imes 1 μ m in size). Scale bars are 10 μ m. (d) Raman spectra in a var2-1 leaf cells were obtained after 5 to 8 h treatment with 100 μM 1, 2, or 3. Three spectra were obtained at different positions in each cite.

and 2. On the other hand, the Raman peak of the diphenylacetylene moiety was not detected in the cell treated with 3 (Fig. 2c and d). Similar results were obtained with root cells (Fig. S10, ESI†). These results indicate that the pLys and TEG components can be applied in planta.

For Raman probes 1 and 2, which showed internalization into plant cells, we carried out a cell viability assay in BY-2 cells (Fig. S11, ESI†). Notably, the cytotoxicity of 1 was higher than that of 2 at the concentration used in this experiment (100 μ M), and TEG-modified probe 2 showed negligible cytotoxicity, which is associated with the different interactions of 1 and 2 with biological membranes. Then, we further investigated the intracellular uptake mechanism of the probe using BY-2 cells in which the intracellular structure is easier to observe. The membrane was stained with an endocytic tracer FM4-64, and the cells were observed by confocal laser scanning microscopy




CSLM imaging of BY-2 cells. (a) Cells treated with each Raman probe were stained with FM4-64 (yellow). (b) Expanded CLSM images for the cells treated with 1, 2, or 3. Arrowheads show vesicles with a size less than 500 nm induced by internalization of 1. Scale bars are 10 µm.

(CLSM) (Fig. 3a). The biological membrane in BY-2 cells was entirely stained in all the cells treated with three probes. In particular, only the cells treated with 1 showed characteristic vesicle structures with a size less than 500 nm stained with FM4-64 in the cytoplasmic area in addition to the other biological membrane (Fig. 3b). We assume that these vesicles are endosomes formed by endocytosis. This result suggests that both probes 1 and 2 are taken into the cells, but their internalization pathways are different.

To investigate whether intracellular internalization of 1 is dependent on endocytosis, we performed two kinds of endocytosis inhibition assays. First, we performed low-temperature treatment to stop the major endocytosis pathway in the cells.²⁶ BY-2 cells were treated with 1 or 2 and incubated at 4 °C for 5 h. Then, the probe signal in the cytoplasm was detected by Raman microscopy. As a result, the probe signal was detected in the cells treated with 2 but not with 1 (Fig. 4a and Fig. S12, ESI†).

Subsequently, the same experiment was conducted by treatment with Wortmannin, an endocytosis inhibitor. Wortmannin is a specific inhibitor of the phosphatidylinositol (PI) 3-kinase that has been shown to block endocytosis in mammalian cells.²⁷ It also has been shown to be active in plant cells, where it inhibits protein sorting to the vacuole through action on both the PI 3- and PI 4-kinases.²⁸ In the cells treated with Wortmannin, the signal of 2 was detected, whereas the signal of 1 was faintly observed (Fig. 4b and Fig. S13, ESI†). Compared with untreated cells, both low-temperature and Wortmannin



Communication

Fig. 4 Raman spectra of BY-2 cells in the cytoplasm. (a) The cells were treated with 1 or 2 at 4 °C. (b) The cells were treated with 1 or 2 in the presence of Wortmannin. The spectra were obtained after 5 h of treatment with 100 μ M **1** or **2**. The cytoplasm area in BY-2 cells (1 μ m \times 1 μ m in size) used for obtaining the averaged Raman spectra are presented by squares in the images. Scale bars are 10 μm.

treatments provided a significant decrease in the Raman intensity at 2215 cm⁻¹ in the BY-2 cells treated with 1 (Fig. S14, ESI†). From these results, it was clear that endocytosis is a dominant pathway for the internalization of 1 into BY-2 cells. In contrast, because 2 was efficiently internalized into the cells even under endocytosis-inhibited conditions, it is assumed that the internalization of 2 is most likely via an energy-independent pathway such as direct penetration.

Hitherto, transmission electron microscopy and atomic force microscopy analyzes have shown that gold nanoparticles coated with pLys or PEG are encapsulated in vesicles and taken up to cells.^{29,30} Nevertheless, the probe containing the PEG component was not applicable to plant cells. It is considered that the PEG-probe caused aggregate formation, so that it could not pass through the cell wall meshwork structure.

Finally, to investigate the intracellular localization of these probes in detail, we performed SRS imaging analysis because of its rapid visualization for the distribution of the Raman probes in cells with high sensitivity. 15,16 We observed BY-2 cells treated with the Raman probes by SRS microscopy. SRS images were obtained at three Raman shifts: 2190, 2235, and 2215 cm⁻¹. The peak intensity at 2215 cm⁻¹ was used to visualize the cells after subtracting the background of the average between 2190 and 2235 cm⁻¹. The obtained images revealed that the three types of Raman probes show different intracellular internalization (Fig. 5).

In the cells treated with 1, a strong signal was detected at a specific site in the cells by SRS imaging. Since we assume that the internalization of 1 proceeds via endocytosis, it is considered that the accumulation site was captured by SRS imaging. Meanwhile, the signal of 2 was detected throughout the cytoplasm. Unlike 1, 2 showed efficient intracellular uptake by direct penetration. The energy-independent direct penetration




Fig. 5 SRS imaging of BY-2 cells treated with Raman probes. BY-2 cells were incubated for 8 h in solution containing Raman probes 1-3 (100 μ M). The images at 2215 cm⁻¹ (left) show the diphenylacetylene moiety of the Raman probe. Background (right) is the images using average intensities of 2190 and 2235 cm⁻¹. The arrowhead shows the putative endosome observed in the cell treated with 1. Scale bars are 10 μm.

pathway probably allows easier accumulation in cells than the energy-dependent endosomal pathway.31 In addition, probe 3 was not detected in any of the cells, which is consistent with the results of Raman microscopy and CLSM shown in Fig. 2 and 3. The result revealed that the structural difference of the polymeric components (pLys or TEG) resulted in the different internalization and localization which play an important role in delivery of exogenous materials into plant cells.

In summary, we demonstrated that the novel alkyne compounds consisting of a diphenylacetylene moiety with a high signal intensity and hydrophilic polymeric components, pLys and TEG. Interestingly, these Raman probes exhibit different intracellular internalization mechanisms depending on the added macromolecular components. The pLys component is considered to interact with biological membranes and dominantly induce endocytosis. In contrast, since TEG shows little interaction with membranes, the probe with the TEG component was proven to be most likely internalized into the cells via direct penetration with negligible cytotoxicity. Based on the findings obtained from this study, practical Raman probes showing high Raman sensitivity, high permeability in plant cells,

RSC Chemical Biology Communication

and low cytotoxicity can be developed by combining diphenylacetylene and specific polymeric components. As an ongoing study, by attaching a reactive functional group to diphenylacetylene, such as a maleimide group, it is possible to bind target cargos, including plasmid DNA, enzyme proteins, and other physiologically active molecules to be delivered intracellularly. The modification with a functional group to the diphenylacetylene moiety is expected to enable tracking of the target cargos by the Raman probe with maintaining the function of pLys or TEG as a membrane-permeable carrier. Visualization of molecules using such Raman probes in plant cells will allow us to obtain deep insights for biological research on the internalization process into plant cells during material delivery.

Conflicts of interest

The authors have no conflicts to declare.

Acknowledgements

We thank Dr W. Sakamoto (Okayama University) for providing the var2-1 mutant plant. This work was supported by IST ERATO Grant Number JPMJER1602, JST CREST Grant Number JPMJCR1872, and JSPS KAKENHI Grant Number JP18K18847.

Notes and references

- 1 B. B. Buchanan, W. Gruissem and R. L. Jones, *Biochemistry* & molecular biology of plants, American Society of Plant Physiologists, Rochville, MD, 2000, pp. 1250-1318.
- 2 A. G. Atanasov, B. Waltenberger, E.-M. Pferschy-Wenzig, T. Linder, C. Wawrosch, P. Uhrin, V. Temml, L. Wang, S. Schwaiger, E. H. Heiss, J. M. Rollinger, D. Schuster, J. M. Breuss, V. Bochkov, M. D. Mihovilovic, B. Kopp, R. Bauer, V. M. Dirsch and H. Stuppner, Biotechnol. Adv., 2015, 33, 1582-1614.
- 3 T. Miyamoto, K. Tsuchiya and K. Numata, Biomacromolecules, 2019, 20, 653-661.
- 4 T. Yoshizumi, K. Oikawa, J.-A. Chuah, Y. Kodama and K. Numata, Biomacromolecules, 2018, 19, 1582-1591.
- 5 M. Fricker, J. Runions and I. Moore, Annu. Rev. Plant Biol., 2006, 57, 79-107.
- 6 A. J. Melveger and R. H. Baughman, J. Polym. Sci., Polym. Phys. Ed., 1973, 11, 603-619.
- 7 H. Yamakoshi, K. Dodo, M. Okada, J. Ando, A. Palonpon, K. Fujita, S. Kawata and M. Sodeoka, J. Am. Chem. Soc., 2011, 133, 6102-6105.
- 8 H. Yamakoshi, K. Dodo, A. Palonpon, J. Ando, K. Fujita, S. Kawata and M. Sodeoka, J. Am. Chem. Soc., 2012, 134, 20681-20689.

- 9 M. Ueda, S. Egoshi, K. Dodo, Y. Ishimaru, H. Yamakoshi, T. Nakano, Y. Takaoka, S. Tsukiji and M. Sodeoka, ACS Cent. Sci., 2017, 3, 462-472.
- 10 G. R. Littlejohn, J. C. Mansfield, D. Parker, R. Lind, S. Perfect, M. Seymour, N. Smirnoff, J. Love and J. Moger, Plant Physiol., 2015, 168, 18-28.
- 11 N. Gierlinger and M. Schwanninger, Spectroscopy, 2007, 21, 69-89.
- 12 N. Gierlinger, T. Keplinger and M. Harrington, Nat. Protoc., 2012, 7, 1694-1708.
- 13 D. R. Chisholm, R. Lamb, T. Pallett, V. Affleck, C. Holden, J. Marrison, P. O'Toole, P. D. Ashton, K. Newling, A. Steffen, A. K. Nelson, C. Mahler, R. Valentine, T. S. Blacker, A. J. Bain, J. Girkin, T. B. Marder, A. Whiting and C. A. Ambler, Chem. Sci., 2019, 10, 4673-4683.
- 14 Y. Chen, J.-Q. Ren, X.-G. Zhang, D.-Y. Wu, A.-G. Shen and J.-M. Hu, Anal. Chem., 2016, 88, 6115-6119.
- 15 L. Wei, F. Hu, Y. Shen, Z. Chen, Y. Yu, C.-C. Lin, M. C. Wang and W. Min, Nat. Methods, 2014, 11, 410-412.
- 16 Y. Ozeki, W. Umemura, Y. Otsuka, S. Satoh, H. Hashimoto, K. Sumimura, N. Nishizawa, K. Fukui and K. Itoh, Nat. Photonics, 2012, 6, 845-851.
- 17 K. Numata, Y. Horii, K. Oikawa, Y. Miyagi, T. Demura and M. Ohtani, Sci. Rep., 2018, 8, 10966.
- 18 A. Salic and T. J. Mitchison, Proc. Natl. Acad. Sci. U. S. A., 2008, 105, 2415-2420.
- 19 T. Nagata, K. Sakamoto and T. Shimizu, In Vitro Cell. Dev. Biol., 2004, 40, 163-166.
- 20 H. Otsuka, Y. Nagasaki and K. Kataoka, Adv. Drug Delivery Rev., 2003, 55, 403-419.
- 21 J. M. Martínez-Zapater, J. Hered., 1993, 84, 138–140.
- 22 M. Chen, J. Hered., 1999, 90, 207-214.
- 23 Y. Kato, E. Miura, K. Ido, K. Ifuku and W. Sakamoto, Plant Physiol., 2009, 151, 1790-1801.
- 24 M. Chen, Y. Choi, D. F. Voytas and S. Rodermel, Plant J., 2000, 22, 303-313.
- 25 W. Sakamoto, A. Zaltsman, Z. Adam and Y. Takahashi, Plant Cell, 2003, 15, 2843-2855.
- 26 T. Wileman, C. Harding and P. Stahl, Biochem. J., 1985, 232, 1-14.
- 27 A. Arcaro and M. P. Wymann, Biochem. J., 1993, 296, 297-301.
- 28 K. Matsuoka, D. C. Bassham, N. V. Raikhel and K. Nakamura, J. Cell Biol., 1995, 130, 1307-1318.
- 29 P. Nativo, I. A. Prior and M. Brust, ACS Nano, 2008, 2, 1639-1644.
- 30 J. K. Vasir and V. Labhasetwar, Biomaterials, 2008, 29,
- 31 S. Trabulo, A. L. Cardoso, M. Mano and M. C. P. De Lima, Pharmaceuticals, 2010, 3, 961-993.